




Fake rate vs normalized χ^2



Duplicates Rate vs normalized χ^2



Pileup rate vs normalized χ^2



Fake rate vs. χ^2

Legend: BQM mkFit5 (blue), BQM mkFit6 (orange), BQM mkFit7 (red), BQM mkFit8 (green), BQM mkFit9 (purple), BQM mkFit10 (brown), BQM mkFit11 (pink), BQM mkFit12 (gray), BQM mkFit13 (light blue), BQM mkFit14 (light orange), BQM mkFit15 (light green), BQM mkFit16 (light purple), BQM mkFit17 (light brown), BQM mkFit18 (light pink), BQM mkFit19 (light gray), BQM mkFit20 (light light blue), BQM mkFit21 (light light orange), BQM mkFit22 (light light green), BQM mkFit23 (light light purple), BQM mkFit24 (light light brown), BQM mkFit25 (light light pink), BQM mkFit26 (light light gray), BQM mkFit27 (light light light blue), BQM mkFit28 (light light light orange), BQM mkFit29 (light light light green), BQM mkFit30 (light light light purple), BQM mkFit31 (light light light brown), BQM mkFit32 (light light light pink), BQM mkFit33 (light light light gray), BQM mkFit34 (light light light light blue), BQM mkFit35 (light light light light orange), BQM mkFit36 (light light light light green), BQM mkFit37 (light light light light purple), BQM mkFit38 (light light light light brown), BQM mkFit39 (light light light light pink), BQM mkFit40 (light light light light gray).

